

Compressor Nebulizer

Model: ME201



MESURE is a high-performance compressor nebulizer which is easy to use with low noise level. It is suitable for continuous use and has two thousand plus hours life service.



ACCESSORIES



AIR TUBE



MOUTH PIECE



ADULT & CHILD
MASKS



AIR FILTER



CUP

SPECIFICATIONS:

Medication Capacity: 6ml

Particle Size: Respirable Fraction 0.5 to 5µm 0.5µm 80%

MMAD: 2.0µm

Average Nebulization Rate: Above 0.2ml/min.

Noise Level: Below 60 dBA

Compressor Pressure Range: 30-45 Psi

Operation Pressure Range: 12-14,5Psi (0.85-1.1 bar)

Liter Flow Range: 8-10 lpm