

**Solutions that flow smoothly!**



**UROLOGY**

# Silicone Foley Catheters - Sterile

Disposable  
100% Silicone  
2way / 3way



## Features:

- Smooth tapered rounded tip.
- Hard valve for easy inflation and deflation of balloon.
- Large smooth eyes.
- Smooth outer surface.
- Strong symmetrical balloon with user friendly deflation and inflation mechanism.
- Fitted with hard plastic valve for easy connection to syringe.
- Colour code sleeve for easy identification.
- Reinforced tip with good catheter flexibility.

| Product code | Products Name                          | Size (Fr/Ch) | Balloon Capacity |
|--------------|----------------------------------------|--------------|------------------|
| OPFCS2106    | Silicone 2-way Silicone Foley Catheter | 6 Fr         | 3 ml             |
| OPFCS2108    | " "                                    | 8 Fr         | 3-5 ml           |
| OPFCS2110    | " "                                    | 10 Fr        | 3-5 ml           |
| OPFCS2112    | " "                                    | 12 Fr        | 5-10ml           |
| OPFCS2114    | " "                                    | 14 Fr        | 5-15ml           |
| OPFCS2116    | " "                                    | 16 Fr        | 5-15ml           |
| OPFCS2118    | " "                                    | 18 Fr        | 5-15ml           |
| OPFCS2120    | " "                                    | 20 Fr        | 30 ml            |
| OPFCS2122    | " "                                    | 22 Fr        | 30 ml            |
| OPFCS2124    | " "                                    | 24 Fr        | 30 ml            |
| OPFCS2126    | " "                                    | 26 Fr        | 30 ml            |
| OPFCS3114    | Silicone 3-way Silicone Foley Catheter | 14 Fr        | 5-15ml           |
| OPFCS3116    | " "                                    | 16 Fr        | 5-15ml           |
| OPFCS3118    | " "                                    | 18 Fr        | 5-15ml           |
| OPFCS3120    | " "                                    | 20 Fr        | 30 ml            |
| OPFCS3122    | " "                                    | 22 Fr        | 30 ml            |
| OPFCS3124    | " "                                    | 24 Fr        | 30 ml            |



STERILE EO

