




"Comfortable Care, Every breath"

NEBULIZER  
MASKS

# Nebulizer Masks



## Features:

- Atomization rate is about 0.4ml/min.
- Drive gas flow is about 4 to 8 L/min.
- Atomization particle  $\leq 5\mu\text{m}$
- Easy-seal, threaded cap and 6cc/8cc capacity jar.
- Anti-spill design prevents loss of medication in any position.
- Jet stays in place unless intentionally removed.
- The product can be transparent green and transparent white.
- The star lumen tubing can ensure oxygen flow even if the tube is kinked, different length of the tubing is available.



Smooth Oxygen Flow

Comfortable Fit

Easy Seal

Product code	Product Specification	Packing
OPNMAE111	Nebulizer Mask with 2.1M tubing Elongated Adult	100pcs/box
OPNMPE112	Nebulizer Mask with 2.1M tubing Elongated Pediatric	100pcs/box
OPNMAS211	Nebulizer Mask with 2.1M tubing Standard Adult	100pcs/box
OPNMPS212	Nebulizer Mask with 2.1M tubing Standard Pediatric	100pcs/box

