




Yankauer and Suction Tubing

Features:

- Made of medical grade PVC, DEHP free available
- Customized with different diameters and lengths
- With different types of connectors
- Step connector of yankauer handle secures strong connectiong with suction connecting tube
- Paper-film pouch packing / blister packing



| OPPY Cat# | Product Description | Pcs/Box | Carton pcs/CTN |
|-----------|---|---------|----------------|
| OPSRYS202 | Suction Catheter Rigid Yankauer w/ Suction Control Valve 1/4 X 3/8IN suction set with 2 mtr tube and handle. ID: 6mm; OD:10mm; Tubing length: 200±5cm; 1 yankauer handle, Flat tip with vent. | 20 | 300 |
| OPSRYS204 | Suction Catheter Rigid Yankauer w/ Suction Control Valve 1/4 X 3/8IN suction set with 3 mtr tube and handle. ID: 7mm; OD:10mm; Tubing length: 300±5cm; 1 yankauer handle, Flat tip with vent. | 20 | 300 |
| OPYSHS102 | Yankauer Handle Plain Tip with Vent | 20 | 300 |



OppyMed LLC 2701 Centerville Road, New Castle County, Willington, Delaware 19808 - USA **E:** info@oppymed.com

EMEA: **OppyMed LLC-FZ** P.O. Box: 9305, 6th Floor, Business Center,
Meydan Hotel Grandstand, Meydan Road, Nad Al Sheba, Dubai, UAE **T:** +971 451 63600 Ext:2313

www.oppymed.com