



**Safe and Secure
Eye Protection
Solution
for Neonatal
Phototherapy**



**NEONATAL
PHOTOTHERAPY
EYEMASK**

NEONATAL PHOTOTHERAPY EYEMASK



OppyMed Phototherapy Eyemask is designed for neonates undergoing phototherapy, offering comfort and eye protection. Its wraparound design secures at two points for a snug fit and prevents movement. The ergonomic eye pad reduces pressure and enhances light protection, while the headband allows blue light to reach the head for continued treatment.



* Ergonomic design for secure, comfortable fit ensure safe and comfortable eye coverage.



* Velcro + elastic fabric, adjustable and easy to use.



* With a soft, non-woven material, absorbent and breathable, reducing friction



* Black shading material, effectively blocks blue light.

Pure cotton non-woven fabric.

Breathable & Absorbent

Reducing friction.

Velcro + Elastic fabric.

Adjustable & Easy to use.

INDICATIONS:

Eye protection for neonatal phototherapy to avoid blue-ray.

TIPS:

- Please consult a Doctor or a Nurse before using this product.
- Please select appropriate head circumference and the size of the protective eye mask.
- Take precautionary measures to patients with allergies and medical eye conditions.

Cat#	Product Description
OPPTES01	Phototherapy Eye shield - Small
OPPTES02	Phototherapy Eye shield - Medium
OPPTES03	Phototherapy Eye shield - Large



OppyMed LLC 2701 Centerville Road, New Castle County,
Wilmington, Delaware 19808 - USA E: info@oppymed.com

www.oppymed.com