





Disposable Warming Blanket

OppyMed Disposable Warming Blankets are designed to maintain patient body temperature before, during, and after procedures, helping reduce the risk of hypothermia. Suitable for use in the waiting room, ICU, operating theatre, and recovery room, they feature flexible, easy-to-use designs, a uniform perforation pattern for consistent convective warming, and are made from durable, soft, radiolucent, latex-free material for enhanced patient comfort and safety.

	Ref #/ Product Description	Size	Pack Size Pc/Box
 Full Body-Adult	OPPWB-FB-A Patient Warming Blanket Full Body-Adult	200 cm x 110 cm (79 in x 43 in)	20 Pc
 Full Upper Body-Adult	OPPWB-UP-B Patient Warming Blanket- Upper Body	215 cm x 100 cm (84 in x 39 in)	40 Pc
 Full Body Pediatric	OPPWB-FB-P Patient Warming Blanket Full Body-Pediatric	150 cm x 100 cm (59 in x 39 in)	30 Pc
 Under Body	OPPWB-UN-B Patient Warming Blanket Under Body	205 cm x 100 cm (81 in x 39 in)	30 Pc
 Under Body Pediatric-Small	OPPWB-UNB-PS Patient Warming Blanket Under Body Pediatric-Small	129 cm x 100 cm (51 in x 39 in)	50 Pc
 Lower Body-Adult	OPPWB-LB-A Patient Warming Blanket Lower Body -Adult	128 cm x 100 cm (50 in x 39 in)	30 Pc